



## Hardness (Total & Calcium) Test Kit, Model HAC-DT

**Product #:**

**2063900**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10 - 4000 mg/L. Contains approximately 100 of each (Total and Calcium) test.

**Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration**

**Rugged titration device is ideal for field use, with no glass burets required**

---

### Specifications

Case Style:	F
Method:	Digital Titrator/EDTA
Number of tests:	100 Ca
	+100 T
Parameter:	Hardness, total & Calcium - As $\text{CaCO}_3$
Platform:	Digital Titrator
Range:	10 - 4,000 mg/L $\text{CaCO}_3$
Smallest increment:	0.1-10
Test Kit Model:	HAC-DT
What's included?:	Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.

---

### What's included?

Includes 6 reagents, digital titrator, delivery tubes, graduated cylinder, flask, manual and carrying case.