



Test kit, hypochlorite, model CN-HRDT

Product #:

2687100

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Hypochlorite Kit, Model CN-HRDT with Digital Titrator

Utilizing the Hach Digital Titrator, this kit is used to determine hypochlorite, high range (Bleach) solution concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 5 - 15% as Cl_2 : 50 - 150 g/L as Cl_2 . Contains 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Designed for on-site analysis of bulk hypochlorite deliveries or bulk storage solutions at drinking water or chloramination facilities

Verify concentration of hypochlorite solutions at time of delivery

Kit contains all required reagents, solutions, and glassware

Specifications

Case Style:	G
Method:	10100
Number of tests:	100
Parameter:	Hypochlorite, high range (Bleach) - As Cl_2
Platform:	Digital Titrator
Range:	5 - 15 % = 50-150 g/L Cl_2
Smallest increment:	0.05%
Test Kit Model:	CN-HRDT