



Test kit, chloride, model CD-DT

Product #: 2063500
ZAR Price: Contact Hach
Ships within 3-5 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Contains mercury. Dispose of according to local, state and federal regulations.

The single-parameter Model CD-DT uses the Hach Digital Titrator with Mercuric Nitrate Titrant Cartridges to determine chloride concentrations. Kit includes Digital Titrator, reagents and labware, manual and carrying case. Range: 10-8,000 mg/L. Contains approximately 100 tests.

Hach Digital Titrator offers higher levels of precision and accuracy than drop count procedures, with reproducibility comparable to a buret titration

Rugged titration device is ideal for field use, with no glass burets required

Specifications

Case Style:	F
Method:	Digital Titrator/ Mercuric nitrate
Number of tests:	100
Parameter:	Chloride - As Cl^- , 10-8,000 mg/L Cl^-
Platform:	Digital Titrator
Range:	10 - 8,000 mg/L Cl^-
Smallest increment:	0.1-20
Test Kit Model:	CD-DT