



## Testkit, glycol visual, model EG-1

**Product #:**

**2186400**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Semiquantitative determination of glycol in water or oil. Visual appearance of colour indicates glycol presence.

Ethylene glycol (anti-freeze) contamination of automotive lubricants is an indication of a cracked engine block. Glycol presence also can indicate leaks in cooling systems. Test Kit Model EG-1 presents a purple color if ethylene glycol and other 1,2 glycols are present in oil or water. The kit contains approximately 25 tests and comes with a 14.4" x 8" x 7" Blue Polypropylene Carrying Case.

**Kit contains all required reagents and apparatus in a rugged carrying case for analysis in the field**

---

### Specifications

Case Style:	G
Method:	Visual/  Appearance of color
Number of tests:	25
Parameter:	Glycol - Ethylene Glycol in Oil or Water
Range:	Presence or absence of glycol, 0- greater than 150 mg/L
Test Kit Model:	EG-1