



## Test kit, sulfate, model SF-1

**Product #:** 225100  
**ZAR Price:** Contact Hach  
**Ships within 3 days**

Widely found in natural waters, sulfate may be at high levels in mine drainage. High concentrations of magnesium or sodium sulfate in water act as cathartics. This test kit provides economical determination of sulfate levels using the Extinction/Turbidimetric methods.

### Economical test kit for determination of sulfate levels using the Extinction/Turbidimetric Method

---

#### Specifications

Case Style:	D
Method:	Extinction/Turbidimetric
Number of tests:	100
Parameter:	Sulfate
Range:	50 - 200 mg/L
Smallest increment:	50
Test Kit Model:	SF-1