



## Test kit, RapidSilver™ visual

Product #: 2674500

**OBSOLETE ITEM**  
This item is no longer available.

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Test Kit Silver, Model RapidSilver. Range 0-50 ppb. Silver As Ag Visual Comparison.

**Provides fast, accurate information on silver discharges**

**Used by surface finishers, photographic film manufacturers and processors, and others to prevent pollution and economic loss**

**Kit directions also include digestion and dilution procedures needed for testing silver strike bath samples**

---

### Specifications

Case Style:	F
Method:	Visual Comparison
Model:	NA
Number of tests:	100
Parameter:	Silver
Range:	5 - 50 ppb Ag
Smallest increment:	5 µg/L
Standard range:	5 - 50 ppb Ag
Test Kit Model:	RapidSilver™