



## Reagent and apparatus set, water in oil

Product #: 2237300

 **Hazardous**

**OBSOLETE ITEM**  
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

### Model WO-1 Test Kit

The Model WO-1 test kit provides critical data in oil quality determinations, oil line operation, and diesel engine maintenance. The kit has proven particularly useful on board ships. Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated.

### Measures low levels of water in oil

Provides critical data in oil quality determinations, oil line operation, and diesel engine maintenance

Reagent is packaged in unit-dose vials, and parts in contact with the sample are disposable, so clean-up is eliminated

---

## Specifications

Case Style:	no case
Method:	Volumetric displacement/  Calcium hydride
Number of tests:	25
Parameter:	Water-in-Oil - As %
Range:	0.05 - 1 %
Range 2:	0 - 10 %
Smallest increment:	0.05%  0.20%
Test Kit Model:	WO-1