



Arsenic Low Range Test Kit

Product #: 2800000
ZAR Price: Contact Hach
Available

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Contains mercury. Dispose of according to local, state and federal regulations.

Hach arsenic test strips provide results within a range of 0 - 500 ppb. This reliable method recovers 100% of inorganic arsenic in samples and allows detection of arsenic levels in the field to assess compliance with the EPA Arsenic Rule.

No lab analysis needed—detect arsenic levels in the field to assess compliance with the EPA Arsenic Rule

Reliable test strip method recovers 100% of inorganic arsenic in samples

Specifications

Case Style:	D
Number of tests:	100
Parameter:	Arsenic, Low Range
Platform:	Test Strip
Range:	0 - 500 ppb
Smallest Increments Steps:	Steps: 0, 10, 30, 50, 70, 300 , 500 ppb
Test Kit Model:	Test Strip†
What's included?:	Includes five reagents, test strips, reaction vessel with cap, measuring scoop (color may vary), instructions and carrying case.

What's included?

Includes five reagents, test strips, reaction vessel with cap, measuring scoop (color may vary), instructions and carrying case.