



EZ Arsenic High Range Test Kit, 0 - 4000 ppb

Product #: 2822800

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Contains mercury. Dispose of according to local, state and federal regulations.

Hach's most economical test kit for determination of arsenic levels.

Simple, economical kit for measuring arsenic levels in the field

Provides fast, reliable results in 20 minutes

Specifications

Case Style:	D
Footnote:	*mg/L unless otherwise noted; ppb = $\mu\text{g/L}$; ppm = mg/L.; gpg = grains per gallon; 1 gpg = 17.1 mg/L or 17.1 ppm.
Number of tests:	100
Parameter:	Arsenic, High Range
Platform:	Test Strip
Range:	0 - 500 ppb As
Range 2:	0 - 4000 ppb As
Smallest increment:	Steps: 0, 10, 25, 50, 250, 500 ppb
Smallest Increments Steps:	Steps: 0; 35; 75; 175; 500; 1,500; 4,000 ppb
Test Kit Model:	Test Strip (EZ)‡
Weight:	1.13 kg