



Be Right™



Quantab Chloride test strips, 300-6000 mg/L, 0.05-1.0 % NaCl , 40 pcs

Product #:

2751340

ZAR Price:

Contact Hach

Ships within 1 week

Ideal for determining quantitative results quickly, simply dip the strip in the sample and compare the length of the colour change to a printed scale, Quantab titrators

Specifications

Method:	Titration Strip
Number of tests:	40
Parameter:	Chloride
Platform:	Test Strip
Range:	300 - 6000 ppm Cl ⁻
Smallest Increments Steps:	100-200 ppm increments
Units:	40 tests