



Test strips, Nitrate 0 - 50 mg/L and Nitrite 0 - 3 mg/L, 25 tests

Product #: 2745425
ZAR Price: Contact Hach
Ships within 1 week

Ranges

Nitrate NO₃-N: 0 - 50 mg/L

Nitrite NO₂-N: 0 - 3 mg/L

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Specifications

Parameter:	Nitrate and Nitrite
Platform:	Test Strip
Range:	Nitrate 0 - 50 mg/L NO ₃ -N mg/L
Range 2:	Nitrite 0 - 3.0 mg/L NO ₂ -N
Smallest increment:	Nitrite: 0, 0.15, 0.3, 1, 1.5, 3 mg/L
Smallest Increments Steps:	Nitrate: 0, 1, 2, 5, 10, 20, 50 mg/L
Test Kit Model:	Test Strips
Units:	25 tests