



Phosphate test strips, 0 - 50 mg/L, 50 tests

Product #: 2757150
ZAR Price: Contact Hach
Ships within 1 week

Easy to use, disposable, and inexpensive

A great way to obtain quick, quantitative answers in the field or in the lab

Specifications

Parameter:	Phosphorus, Orthophosphate (reactive) - as PO_4
Platform:	Test Strip
Range:	0 - 50 ppm PO_4
Smallest Increments Steps:	Steps: 0, 5, 15, 30, 50
Test Kit Model:	Test Strips
Units:	50 tests