



Filters, Membrane, 47-mm, 0.45-µm, gridded, sterile, Millipore

Product #:

2936100

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Sterile membrane filters for use in microbiology membrane filtration methods.

Specifications

Diameter:	47 mm mm
Material:	Mixed-Cellulose Ester
Pore Size:	0.45 µm
Quantity:	150