



Be Right™



Bottle, Sampling, Polycarbonate, 1000 mL, 4/pk

| | |
|-------------------|----------------|
| Product #: | 2324663 |
| ZAR Price: | Contact Hach |
| Available | |

Transparent polycarbonate bottle with cap features wide mouth (3.8 cm/1.48" diameter) opening for easy filling and square design for maximum space efficiency. Autoclavable with molded-in graduations (do not engage threads for autoclave cycle.)

Specifications

| | |
|----------------|---------------|
| Body Material: | Polycarbonate |
| Capacity: | 1000 mL |
| Quantity: | 4 /pk |