



## Filter, acrocap, pore size 0.2 µm

<b>Product #:</b>	<b>2323810</b>
ZAR Price:	Contact Hach
Available	

For sterilization or clarification of up to three liters of culture media or other aqueous solutions. Ideal for sterilization of solutions that tend to foam when filtered under vacuum. Inlet accepts 6.4 mm (0.25") tubing or most male luer fittings.

---

### Specifications

Diameter:	NA
Filter Type:	Filter, acrocap
Packing:	NA
Pore Size:	0.2 µm