



## Filter, acrodisc, pore size 0.2 µm

<b>Product #:</b>	<b>2415600</b>
ZAR Price:	Contact Hach
Available	

For filtration of dilute biological fluids. Sterilized by gamma irradiation to eliminate potential contamination by EtO residuals. Inlet/Outlet connections: 13/25 mm: Female Luer-Lok inlet, standard male slip luer outlet. Pore size: 0.2 µm. Diameter: 25mm.

---

### Specifications

Diameter:	NA
Filter Type:	Filter, acrodisc
Material:	Polysulfone
Packing:	NA
Pore Size:	0.2 µm