



Microfil funnel & filter, 100 mL, 0.45 µm 47 mm white gridded (millipore)

Product #:	2586300
ZAR Price:	Contact Hach
Ships within 1 week	

Designed to eliminate cleaning and sanitization of filter funnels between uses. Graduated (20, 50, and 100 mL) 100 mL funnels are sterile and ready-to-use. Funnel creates leak-proof seal when seated onto Microfil Stand Alone Support. Excellent for testing water samples from potable, non-potable, source, ambient, industrial, food, beverage, pharmaceutical and wastewater sources. Microfil funnels come with 150 Millipore filters (white; gridded; cellulose and mixed ester composition; 0.45 µm pore size; 47 mm diameter; sterile; individually packed). Microfil V includes integrated filter for added convenience and can be used with standard No. 8 stopper on any filter holder or manifold.

Specifications

Capacity:	100 mL
Diameter:	47 mm
Material:	Mixed-Cellulose Ester
Pore Size:	0.45 µm
Quantity:	150