



Microfil V Filtration Device

Product #: 2508900

OBSOLETE ITEM
This item is no longer available.

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Designed to eliminate cleaning and sanitization of filter funnels between uses. Graduated (20, 50, and 100 mL) 100 mL funnels are sterile and ready-to-use. Funnel creates leak-proof seal when seated onto Microfil Stand Alone Support. Excellent for testing water samples from potable, non-potable, source, ambient, industrial, food, beverage, pharmaceutical and wastewater sources. Microfil funnels come with 72 Millipore filters (white; gridded; cellulose and mixed ester composition; 0.45 µm pore size; 47 mm diameter; sterile; individually packed). Microfil V includes integrated filter for added convenience and can be used with standard No. 8 stopper on any filter holder or manifold.

Specifications

Capacity:	100 mL
Diameter:	47 mm
Material:	Mixed-Cellulose Ester
Pore Size:	0.45 µm
Quantity:	1