



## Filter funnel, micro

**Product #:**

**2619000**

ZAR Price:

Contact Hach

Ships within 6-8 weeks

Feature a completely assembled filtration unit: Funnel, filter, and petri dish. Filter: 47 mm, white with gridded surface, 0.45 µm GN-6 Metrical media Funnel: Polypropylene Petri Dish: polystyrene lid and polypropylene base Fits 28315-00, 28461-00 and 28461-50 manifolds. Sterile and individually bagged.

---

### Specifications

Diameter:	47 mm
Material:	Mixed-Cellulose Ester
Pore Size:	0.45 µm
Quantity:	50