



Membrane filters millipore S-pak type HA, 47 mm, 0.45 um pore size, pk/600

Product #:

2936500

ZAR Price:

Contact Hach

Ships within 1 week

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

For membrane filtration in microbiological analysis of aqueous samples. "HA" signifies certified compliance with Standard Methods criteria with pore size of 0.45 um. White, cellulose and mixed-ester composition. S-Pak (filter only) is sterile. Diameter: 4

Specifications

Diameter:	47 mm mm
Material:	Mixed-Cellulose Ester
Pore Size:	0.45 µm
Quantity:	600