



## Filter membrane, EZ-Pak, 0.22 micron, 47 mm, white, gridded, 600/pk (millipore)

Product #:

2916900

ZAR Price:  
Available

Contact Hach



**Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Dispenser provides superior convenience and assures sterility. Unique one press on the lever and a sterile membrane is presented (horizontally and inverted) for efficient and convenient assembly. The EZ Pak membranes come as a band of membranes that easily loads into the dispenser. Packaging materials are wound into the dispenser ensuring a clear operating area and quick change of bands. The membranes are Type HA, which signifies compliance with Standard Methods criteria.

---

### Specifications

Diameter:	47 mm mm
Material:	Mixed-Cellulose Ester
Pore Size:	0.22 µm
Quantity:	600