



Filter, membrane, pore size 0.2 µm, diameter 47 mm, 200/pk (Gelman)

Product #:

2618800

ZAR Price:

Contact Hach

Ships within 2 weeks

For isolation and enumeration of stressed organisms (*Pseudomonas* species). High retention and flow rate. Certified for use in membrane filtration analysis. White polyethersulfone media. STerile and individually packaged. Pore size: 0.2 µm. Diameter: 47mm.

Specifications

Diameter:	47 mm
Filter Type:	Filter
Material:	Polyethersulfone
Pore Size:	0.2 µm
Quantity:	200 /pk