



Sample vial, 120 mL w/ thiosulfate & custody tie, 200/pk

| | |
|---------------------|----------------|
| Product #: | 2634001 |
| ZAR Price: | Contact Hach |
| Ships within 1 week | |

Sample Vial with Thiosulfate and Custody Tie

120 mL sterile, hinged-cap vials meet all EPA requirements for microbiology sampling. Collection, incubation, and accurate reading can all be done in the vial without a transfer step. 100 mL raised fill line. Polypropylene. (Please Note: These sample vials are not autoclaveable)

Specifications

| | |
|-------------------|---------------------------------|
| Capacity: | fill line at 100 mL (+/-2.5 mL) |
| Diameter: | 43mm |
| Length: | 93mm |
| Quantity: | 200 /pk |
| Special Features: | Sample Vial, with Thiosulfate |
| Unit: | 200 /pk |