



Filter, Gelman, Finished Water

Product #:	2811200
ZAR Price:	Contact Hach
Available	

Use the Envirocheck Filter to concentrate living organisms (Giardia and Cryptosporidium) and non-living particulates. Vinyl caps fit over the hosebarbed inlet and outlet ports, which form a watertight seal. Pore size: 30 psid @ ambient 20-25°C. Dimensions: 6 cm (2.4") diameter, 21 cm (8.5") length. Pall Life Sciences (Gelman).

Specifications

Quantity:	1 /pk
-----------	-------