



## Dish, Petri, 15 x 100 mm

<b>Product #:</b>	<b>2178900</b>
ZAR Price:	Contact Hach
Ships within 1 week	

Widely used in microbiology, these standard-size polystyrene petri dishes are presterilized, stackable, and disposable. Dishes are packaged 24 per sleeve, with a case containing 25 sleeves (600 dishes total).

**Standard-size polystyrene petri dishes**

**Presterilized, stackable, and disposable**

**600 dishes total**

---

### Specifications

Quantity:	600
-----------	-----