



Filter holder, for vacuum filtration, glass

Product #:	234000
ZAR Price:	Contact Hach
Ships within 1 week	

Beaker-style funnel (300 mL) with 47 mm membrane support for filtration of aqueous samples. Fritted borosilicate glass membrane support for even filtration across surface of membrane. Funnel attaches evenly to base via aluminum clamp. Filter flask sold separately.

Specifications

Diameter:	47 mm
Quantity:	1
What's included?:	Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).

What's included?

Borosilicate glass funnel and base; anodized aluminum clamp; support screen (sintered glass); rubber stopper (No. 8).