



Incubator, BOD, Compact, Model 205, 220/240 VAC

Product #: 2616202

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Suggested replacements:
• 2636202

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Hach's compact Model 205 BOD incubator exceeds the incubation criteria for biochemical oxygen demand (BOD) analysis, and easily fits under standard lab benches. Digital temperature display with high-limit warning light. Capacity: 59 standard, 300 mL BOD bottles (3.1 ft³) or one BODTrak. Internal dimensions are 29.5 x 16 x 15.5 inches (H x W x D). External dimensions are 38 x 18.5 x 19 inches (H x W x D). Temperature Range: 5 - 45 °C (Accuracy: ±0.5 °C @ 20 °C)

Specifications

Accuracy:	± 0.5 °C at 20 °C Temperature Uniformity
Capacity:	59 BOD Bottles
Electrical Specifications:	220 Vac, 50/60 Hz
Manual Languages:	English
Method:	8043
Temperature Range:	5 - 45 °C
Weight:	86.2 kg Shipping Weight
	27.67 kg Net Weight
What's included?:	Incubator and 2 shelves

What's included?

Incubator and 2 shelves