



Cooler, Sample

Product #:

2568500

ZAR Price:

Contact Hach

Ships within 3-5 weeks

9 qt. Igloo cooler ideal for transport of temperature sensitive samples. Detachable lid reverses to a tray and cup holder.

Interior Dimensions: (L x W x H) 10.25 x 7.38 x 6.94 inches

Interior Dimensions: (L x W x H) 26 x 18 x 17.5 cm

Exterior Dimensions: (L x W x H) 13 x 9.5 x 8.75 inches

Exterior Dimensions: (L x W x H) 33 x 24 x 22cm

Specifications

Dimensions (H x W x D): 176.276 mm x 187.452 mm x 260.35 mm

Quantity: 1

Weight: 1.32 kg