



Thermometer, 30 to 50 C°, 250 mm

Product #: 2999200
ZAR Price: Contact Hach
Available

Non-mercury thermometer, 30 - 50 degrees C Length: 250 mm Immersion: 76 mm Safety coated.

Blue fill solution on white background provides excellent visibility.

Non-mercury thermometer, 30 - 50 degrees C Length: 250 mm Immersion: 76 mm Safety coated.

Blue fill solution on white background provides excellent visibility.

Specifications

Length:	250 mm
Quantity:	1
Temperature Range:	30 °C - 50 °C
Thermometer Divisions:	0.1 °C
Thermometer Fluid Type:	Non-Mercury, blue