



Digital Colony Counter Pen

Product #:	25045000
ZAR Price:	Contact Hach
Ships within 1 week	

Counter Pen is a combination marker and digital counter that can be used to count cell growth on petri dishes, inventory parts, and tally samples received. To ensure an accurate tally, the counter beeps each time a mark is made, and the easy-to-read LCD displays the total (up to 99,999). Pen features on/off switch and comes with cap.

Count cell growth on petri dishes

Inventory parts

Tally samples received

Specifications

Material:	Durable ABS plastic
Quantity:	1