



Bottle rack, plastic-coated, 15 compartments

Product #: 2551800

 **Hazardous**

OBSOLETE ITEM
This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Plastic-coated wire rack provides corrosion resistant transport system for sample bottles and bags. Racks feature 15 or 20 compartments with double-grid bottom to prevent loss of samples during transport. 15 compartment rack is ideal for Whirl-Pak sample bags (4 oz). 20 compartment rack fits most bottles associated with standard drinking water sampling. Autoclavable.

Specifications

Quantity:	1
Special Features:	Rack, 15 Bottle Compartments