



m-FC Broth, Dehydrated 500 g

Product #: 1462434

OBSOLETE ITEM
This item is no longer available.

m-FC dehydrated Broth Base. 500 g.

Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

Specifications

Media Type:	m-FC Broth
Quantity:	500 g