



## Bottle, sampling, presterilized, round, 50/pk

<b>Product #:</b>	<b>2495050</b>
ZAR Price:	Contact Hach
Ships within 3 days	

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Work well with P-A Broth Ampules and PathoScreen Medium.

---

### Specifications

Dechlorinating Agent?:	No
Quantity:	50/pk