



## Bottle, sampling, presterilized, round, 12/pk

Product #:

2495012

 **Hazardous**

### **OBSOLETE ITEM**

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Polystyrene, disposable, transparent, non-fluorescent, break-resistant, presterilized bottles with caps. Bottles have 100 mL fill-to line; capacity of 175 mL. Bottles used with medium containing MUG or other fluorescing reagent. Works well with P-A Broth Ampules and PathoScreen Medium.

---

## Specifications

Dechlorinating Agent?:

No

Quantity:

12/pk