



Severn Trent Portable Luminometer

Product #:

2887000



Hazardous

OBSOLETE ITEM

This item is no longer available.

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

The Severn Trent Portable Luminometer measures and records relative light units (RLUs) in aqueous samples. The luminometer is made for use under extreme field conditions, and the components are rugged, easy to use, and reliable. The luminometer has stored programs for both Eclox and luminescent bacteria. The luminometer is used in the Hach Eclox Rapid Response Water Kit.

Ruggedized for use in all conditions

The Severn Trent luminometer has ruggedized components and is powered by commonly available (AA) batteries for ultimate reliability

Simple to operate

Stored programs, simple interface, and large readout make the instrument easy to use for a variety of tests

Instrument only - reagents not included

The Severn Trent portable luminometer is used in the Hach Eclox Rapid Response Water Kit, available separately. This product contains the luminometer only.

Measures 0 - 2,000 RLUs

The Severn Trent portable luminometer measures 0 - 2,000 Relative Light Units (RLU)