



MEL/P-A safe drinking water lab

Product #: 2569600
ZAR Price: Contact Hach
Ships within 3-5 weeks

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Contains mercury. Dispose of according to local, state and federal regulations.

MEL/P-A Safe Drinking Water Lab

Economical portable lab monitors key water quality parameters plus Total Coliforms and E. Coli (Presence-Absence Method). MEL P/A Safe Drinking Water Lab includes portable incubator, reagent sets, portable long-wave UV lamp, manual and apparatus.

Consumables only

Choose media and order seperately

Nitrate and Chlorine determined using visual color discs

Portable incubator used to determine total and fecal coliforms directly in the field

Specifications

Sensitivity: CFU/mL