



## PathoScreen Field Test Kit

**Product #:**

**2859100**

**ZAR Price:**

Contact Hach

**Ships within 3 days**

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

PathoScreen Field Test Kit. Enough media for 100 presence-absence tests using 20 mL sample, or 20 MPN tests. Reusable glass vials, and foil-pillowed media, increases shelflife. No refrigeration required for media.

---

### Specifications

Quantity: 1 each

What's included?: Includes Pathoscreen powder pillows, Bleach Solution, Instruction Sheet, and Carrying Case.

---

### What's included?

Includes Pathoscreen powder pillows, Bleach Solution, Instruction Sheet, and Carrying Case.