



## **Hach Modified Colitag Sterile 120 mL Sample Bottles, 100/pk, with Shrink-banded, Polystyrene**

<b>Product #:</b>	<b>8888006</b>
ZAR Price:	Contact Hach
Ships within 2 weeks	

### **Hach Colitag presence-absence tests deliver the fastest EPA approved results in the market**

Hach Modified Colitag Sterile 120 mL Sample Bottles with sodium thiosulfate tablet. Can be used with all Hach Modified Colitag 16-hour presence/absence tests. 100/pk. Hach Modified Colitag offers the following benefits:

#### **The Fastest Presence/Absence Test on the Market at 16 Hours**

Fastest EPA-approved test for coliforms and *E. coli*. Test as late as 5:00 PM and still have your results by 9:00 AM the next day.

#### **One Product that Supports 16- or 22-hour Test Methods**

Allows for a "standard" test that requires less set-up, or a

"rapid" test when necessary. One media for both tests reduces

complexity in certification and inventory.

#### **Results are Valid for 48 Hours**

**(Please Note: When transported in hot weather, the tab may melt. This does not affect product quality or function).**

Saturday's tests can be read on Monday—no need to work all

weekend. Results stay valid for 800% longer than the competition.

---

### **Specifications**

Material:	Plastic
Quantity:	100 /pk