



## Plastic ampules, pseudomonas broth, pk/50

<b>Product #:</b>	<b>2812250</b>
ZAR Price:	Contact Hach
Available	

Pseudomonas broth plastic ampules, 2 mL. For detection of Pseudomonas species. pk/50.

### Complete Solutions for Testing Recreational Water

Hach offers a variety of media and apparatus to meet all of your recreational testing needs.

### Complete Solutions for Testing Wastewater

Hach has prepared media and convenient apparatus to make wastewater and non-potable water testing fast and easy.

### Complete Solutions for Testing Drinking Water

Whichever method a customer chooses for examining drinking water samples, Hach has prepared media and convenient apparatus to make testing easy and fast.

---

## Specifications

Material:	Plastic
Media Type:	Pseudomonas Broth
Quantity:	50 /pk