



## Sample transport kit

Product #:

2568700

### **OBSOLETE ITEM**

This item is no longer available.



Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

Contains everything needed for collecting field samples, keeping them cold during transport. Includes Rubbermaid Cooler, plastic-coated rack with 15 compartments for sampling bags, 100 Whirl-Pak sampling bags (177 mL) with dechlorinating agent, one leak-proof refrigerant pack. Dimensions: 42.4 x 28.4 x 58.6 cm (16.7 x 11.2 x 23.1")