



Formazin Turbidity Standard, 4000 NTU, 500 mL

Product #:	246149
ZAR Price:	Contact Hach
Ships within 3 days	

Turbidity standard can be diluted to be used with any of Hach's turbidimeters.

Specifications

Concentration:	4000 NTU
Measuring range:	N/A
Parameter:	Turbidity
Quantity:	500 mL 500 mL
Storage conditions:	Store in a cool, dark place, away from direct sunlight

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)