



Stablcal Certified calibration standard 800 mNTU for Filtertrak 660

Product #: 2788453
ZAR Price: Contact Hach
Ships within 3 days

For 1-point calibration and 2-point calibration. 1 L bottle.

Specifications

Concentration:	800 mNTU
Measuring range:	N/A
Parameter:	Turbidity
Quantity:	1 L
Storage conditions:	Store between 5 - 25 °C, away from direct sunlight

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)