



## Free Chlorine Reagent Set for chlorine analyser CL17/CL17sc

**Product #:**

**2556900**

**ZAR Price:**  
Available

Contact Hach

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accommodate these items.

### Formulated to meet the strictest quality standards

The Free Chlorine Reagent set is used with the CL17/CL17sc analyser as well as bulk analytical methods. Hach reagents provide a guarantee of accurate, consistent results.

Free Chlorine Reagent Set Consists of:

- One 24-gram bottle of DPD compound
- One 473-mL bottle of free chlorine indicator solution
- One 473-mL bottle of free chlorine buffer solution

\*This product has not been evaluated to test for chlorine and chloramines in medical applications in the United States.

### Quality, Reliability

Hach Free Chlorine Reagent Set (DPD Indicator, Free Chlorine Indicator and Buffer Solutions) is specifically designed and produced to work with Hach CL17 Analyzers.

### Simple, Predictable Maintenance

The CL17 is designed to require replacement of reagents on a predictable schedule under typical use. Maintenance can be performed in 15 minutes, and no special tools are required.

### Reagent Delivery Program

Want to simplify reagent reordering? Ask your Hach representative about our simple and easy automatic delivery program for Hach CL17 reagents and reagent sets.

---

## Specifications

|                     |   |
|---------------------|---|
| Instrument:         | CL17/CL17sc   |
| Method:             | 10059   |
| Number of tests:    | 450   |
| Parameter:          | Chlorine, Free  |
| Platform:           | Solution  |
| Range:              | 0.03 - 5.00 mg/L $\text{Cl}_2$ for CL17, 0.03 - 10.00 mg/L $\text{Cl}_2$ for CL17sc             |
| Storage conditions: | Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture. |

|                   |  |
|-------------------|--|
| What's included?: | DPD compound, 24 g                       |
|                   | Free chlorine indicator solution, 473 mL |
|                   | Free chlorine buffer solution, 473 mL    |

---

**What's included?**

DPD compound, 24 gFree chlorine indicator solution, 473 mLFree chlorine buffer solution, 473 mL

---

**Required Accessories**

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)