



Chlorine Standard, 100-130 mg / L 20 mL Amp

Product #:	2835801
ZAR Price:	Contact Hach
Available	

This chlorine standard can be diluted down to verify the accuracy of the CL17 chlorine analyser or create a separate calibration curve. Under normal circumstances, a separate calibration curve is not required for the CL17. Package of one 20 mL ampule.

Specifications

Instrument:	CL17
Measuring range:	N/A
Storage conditions:	Keep refrigerated.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)