



Electrolyte for 9185sc Sensor, 100 mL

Product #:	09185=A=3600
ZAR Price:	Contact Hach
Ships within 1 week	

Replacement electrolyte for the 9185sc Ozone Sensor.

Specifications

Instrument:	9185sc
Measuring range:	N/A
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)