



Reagent Set, APA6000 Low Range Ammonia/ Monochloramine Analyser

Product #:

6001400

ZAR Price:

Contact Hach

Ships within 3 days



Hazardous

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Specifications

Instrument:	APA6000
Measuring range:	N/A
Parameter:	Ammonia/Monochloramine
Range:	Ammonia 0.02 - 2.0 mg/L NH_3 0.1 - 10.0 mg/L as Cl_2 Monochloramine 0.02 - 2 mg/L as N, 0.1 - 10.0 mg/L as Cl_2
Special Features:	Reagent Set, 4 x 1L
Unit:	4 pk

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)