



APA6000 Ammonia/Monochloramine, Standard 1, 0 mg/L, NH₃, 1 L

Product #: 2776653
ZAR Price: Contact Hach
Available

This standard is used as part of the standard kit for both the APA6000 Ammonia and Ammonia/Monochloramine instruments.

Specifications

Instrument:	APA6000
Measuring range:	N/A
Parameter:	Ammonia/Monochloramine
Storage conditions:	Keep containers tightly closed in a dry, cool and well-ventilated place. Protect from moisture.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)