



## 932100\_Kit-Smart-Probe Adapter-TZ

<b>Product #:</b>	<b>9321000</b>
ZAR Price:	Contact Hach

The Smart Probe Adapter for Hach 5500sc and 961x analyzers allows connection of Hach digital sensors, allowing even greater process control

**5500sc/961x Analyzers can be connected to Hach digital sensors, eliminating the need for additional controller**

**5500sc Ammonia Monochloramine can now be used in conjunction with the CL17sc, CL10sc, or other digital (sc) sensors**