



5500sc Ammonia Monochloramine Reagent 3, 1 L

Product #:

25236000

ZAR Price:

Contact Hach

Ships within 3 days

 **Hazardous**

Items with this mark may be considered hazardous under some shipping conditions.

If necessary, we will change your selected shipping method to accomodate these items.

Reagent 3 for 5500sc AMC Analyzer, 1 L

Specifications

Instrument:	5500sc
Measuring range:	N/A
Parameter:	Ammonia/Monochloramine
Storage conditions:	Keep container tightly closed in a dry and well-ventilated place. Keep in properly labeled containers.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)