



5500sc Ammonia Monochloramine Standard 1 (0 ppm NH₃), 2 L

Product #: 25237000
ZAR Price: Contact Hach
Ships within 3 days

Standard 1 for 5500sc AMC Analyzer, 2 L

Specifications

Instrument:	5500sc
Measuring range:	N/A
Parameter:	Ammonia/Monochloramine
Storage conditions:	Keep container tightly closed in a dry and well-ventilated place. Keep in properly labeled containers.

Required Accessories

- Gloves, Disposable, Powder Free, Nitrile, Medium (Item 2550502)
- Gloves, Disposable, Powder Free, Nitrile, X-Large (Item 2550504)
- Gloves, Disposable, Powder Free, Nitrile, Small (Item 2550501)
- Gloves, Disposable, Powder Free, Nitrile, Large (Item 2550503)